pierre cardin
PARIS

pierre cardin
PARIS

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pierre cardin
PARIS

Safety Precautions

To reduce the risk of malfunction, electric shock, injury, fire, death and damage to equipment or property, please observe the following safety precautions.

Explanation of Symbols

⚠ Danger: will result in serious injury or death of a potential hazard.	⚠ Warning: denotes a potential hazard that could result in serious injury or death.	⚠ Caution: denotes a hazard that could result in minor injury.
---	---	--

❗ This symbol is used to alert users to a specific operating procedure that must be followed in order to operate the unit safely.

⚠ Danger ⚠ Warning

⌛ This symbol is used to alert users to a specific operating procedure that must be followed in order to operate the unit safely.

❗ This symbol is used to alert to a specific operating procedure that must be followed in order to operate the unit safely.

Abnormal or failure cases
The main unit or charger is deformed or abnormally hot.
Burning smell from the toothbrush or charger.
There is abnormal sound during use or charge the toothbrush.
Always ensure the appliance is operated on an electric power source that match the rated voltage. Fully insert the power plug.
Failure may cause fire or electric shock.
The product is not recommended to be played by children neither the disabled person in case of accident or injury.
Replace the charger if the cord or adaptor is damaged.
Remove the battery before disposing of the toothbrush.
Battery must be disposed of carefully.

General Description

Brush Head
On/Off & Function Switch
Fast Middle Slow Rotation
Charging Indicator
Type-C Charge Port

Storing and Charging the Toothbrush

- Charge the device when the toothbrush stops operating; the first time usage over three months without usage or during usage.
- Temperature for charging is 9°C-35°C. The battery may not be charged properly under extreme low or above high temperatures.
- Keep the main unit attached on the charger to maintain the full battery capacity.
- The toothbrush maybe hot during usage or charging.

Charging Tips

- Charging time: 4-6 hours.
- Charging status: indicator flashes.
- Full charged: indicator is normally on.
- Charge temperature: 9°C-35°C.
- Indicator may not light until 6-10 minutes in charging after long time storage.

Maintenance Instructions

Clean the Brush Head
Rinse to remove residues after every brushing and dry it naturally.

Clean the Handle
Rinse to remove residues after every brushing and dry with a soft cloth. Do not immerse the toothbrush into water. Do not clean the toothbrush with water higher than 40°C or corrosive cleanser.

Replace the Brush Head
Professional dentist recommends replacing brush head every 3 months for the best brushing results. Also replace the brush head when the bristles start to be outwards.

Start

Carefully insert the toothbrush head into the output shaft.

Wet the bristles and apply a certain amount of toothpaste.

Put it into your mouth before switching on in case of toothpaste splash. Then press the power button again to select a needed mode.

Guide your toothbrush to move slowly between the teeth, so that each one is cleaned.

Move the bristles slowly in circular motion. Do not scrub teeth, but apply light pressure to let it do the brushing for you.

Gargle your mouth after brushing.

Brushing Time

The timer features 30-second intervals for the full dental recommended 2 minutes of brushing. There is a short pause after each interval reminding to switch to a different area.

Function

2-Min timer
Weak standby time
Inductive charging
Non-slip Handle
IPX7 waterproof
Three working modes

Key Function

ON/OFF Switch 1: Short press control to turn on/off the electric toothbrush.
Function Switch 2: Long press to switch to 3 cleaning modes.

Specification

Item	Rotary electric toothbrush
Power	1,7W DC5V 1A 50-60Hz
Battery	600mAh, 2mins smart timing, 30s reminder
Charging time	4h
Accessory	2 pcs brush-head, 1pc charging cable, 1pc manual, 1pc gift box, 1 travel case
Speed	1st 4500 2nd 6500 3rd 8800

Troubleshooting

Problem	Diagnosis	Solutions
Rotation weakened	1.The bristle is pressed too hard 2.The brush head is fixed too tight	1.Do not press the bristle with too much force 2.Reconnect the brush head with some space in-between
The toothbrush can only run for a short time	1.The toothbrush has been unused for more than 6 months 2.Insufficient charge time	1.Make sure the toothbrush is placed vertically on the charge seat and the charge pin is inserted 2.Charge till the indicator blinks
The toothbrush does not work at all	1. Inner mechanics failure 2.Battery has almost reached its lifetime	Contact your dealer for further inspection

Important Instructions

GENERAL NOTICES

- This appliance is not intended for use by people (including children) with reduced physical sensory or metal capabilities, or lack of experience and knowledge.
- Children between the age of 6-10 must be given supervision or instruction concerning use of the appliance by a person responsible for their safety. Children under the age of 6 are not recommended for use.
- This appliance is only for personal use. It is not intended for use in dental clinics by multiple patients.
- This appliance is for cleaning teeth, gums and tongue only. Do not use it for other purposes.
- Charge indicator may not come to work instantly when you charge after appliance has been unused for a long period of time. It will switch back on in 6-10 minutes.
- It is normal to feel slight heat of the handle during charging process.
- Clean the pins of the power adapter after a certain period of time to avoid overheating and electrocution that may caused by low insulation due to dust accumulation.
- Store the appliance away from direct sunlight where it is well-ventilated and dry. Temperature at 7°C-40°C.
- Do not expose surface of the appliance to certain substance (gel oil, soap, etc.) for a long time to avoid corrosion.
- Disposal of the appliance and its built-in battery should comply with your local regulation.

MEDICAL WARNINGS

It is normal for first-time users to experience slight bleeding and discomfort for their initial week of use. If bleeding or discomfort continues to occur after 3 weeks of use, discontinue and contact your dental professional. If you have had oral or gum surgery in the previous 2 months, implanted dental prosthesis, consult your dental professional prior to use. If you are applying foam or other implanted device, get professional advice with physician in advance. Get related advices if you still have other medical concern.

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Teknikata 4, Diğ. Tic. Ltd. Şti.